



PRECISION BALANCE BBAL-219

PRECISION BALANCE BBAL-219

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

BBAL-219 PRECISION BALANCE



Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing
4 conversion units: oz,ct,lb,g

SPECIFICATIONS

| | |
|-----------------------|--|
| Model | BBAL-219 |
| Capacity (g) | 1100 g |
| Resolution (g) | 0.01 g |
| Min Weighing (g) | 0.04 g |
| Stable Time | ≤2 sec |
| Pan Size | Φ 130 mm |
| Operation Temperature | 5-25°C |
| Repeatability | 0.02 g |
| Liner | 0.03 g |
| Cal. Weight | 500 g |
| Packing Size | Inside: 310x230x130 mm Outside: 660x470x320 mm |
| Gross Weight | 2 / 22 kg |

OPTIONAL ACCESSORIES

| Accessory Code | Name |
|----------------|--------------------------|
| 3600606006 | RS232C & RS485 Interface |
| 3600606007 | Printer |
| 3600715006 | Windbreak mask |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com