



PRECISION BALANCE BBAL-225

PRECISION BALANCE BBAL-225

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

BBAL-225 PRECISION BALANCE



Aluminium shell, Anti-jamming, Anti-radiation

LED digital display

The multiple functions modes with tare, counting, unit conversion and weighing

SPECIFICATIONS

Model	BBAL-225
Capacity (g)	5100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	3 sec
Display	LED / LCD
Sensor	Load Cell
Pan Size	180x180 mm
Operation Temperature	5-25°C
Repeatability	±0.03 g
Cal. Weight	2000 g
Gross Weight	5 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer
3600725006	USB



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com