



## WEIGHING SCALE BBAL-413

## WEIGHING SCALE BBAL-413

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory.

Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research.

Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

## BBAL-413 WEIGHING SCALE



Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing

## SPECIFICATIONS

Model	BBAL-413
Capacity (g)	210 kg
Resolution (g)	10 g
Min Weighing (g)	50 g
Stable Time	≤2 sec
Display	LED
Sensor	Load Cell
Pan Size	410x510 mm / 450x460 mm / 600x800 mm
Operation Temperature	5-25°C
Repeatability	±20 g
Liner	±30 g
Packing Size	420x170x850 mm
Gross Weight	13 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)