



## WEIGHING SCALE BBAL-414

## WEIGHING SCALE BBAL-414

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory.

Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research.

Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

## BBAL-414 WEIGHING SCALE



Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing

## SPECIFICATIONS

Model	BBAL-414
Capacity (g)	310 kg
Resolution (g)	10 g
Min Weighing (g)	50 g
Stable Time	≤2 sec
Display	LED
Sensor	Load Cell
Pan Size	410x510 mm / 450x460 mm / 600x800 mm
Operation Temperature	5-25°C
Repeatability	±20 g
Liner	±30 g
Packing Size	420x170x850 mm
Gross Weight	13 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)