



## SHAKING WATER BATH BBSW-3001

# SHAKING WATER BATH BBSW-3001

## BBSW-3001 SHAKING WATER BATH



A stamping molding stainless steel tank, easy to clean  
LCD screen, multiple data display with timing function, easy to operate  
Stainless steel shelves cover heater and sensor to avoid damage during using  
Once-forming stainless steel lid  
Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time  
Independent temperature-limiting alarm system  
temperature error alarm  
Test tube holder can be placed

## SPECIFICATIONS

Model	BBSW-3001
Capacity	33 L
Temperature range	RT+5-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1250 W
Interior Dimension	438Wx310Dx250H mm
Exterior Dimension	643Wx350Dx353H mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)