



SHAKING WATER BATH BBSW-3005

SHAKING WATER BATH BBSW-3005

BBSW-3005 SHAKING WATER BATH



A stamping molding stainless steel tank, easy to clean
LCD screen, multiple data display with timing function, easy to operate
Stainless steel shelves cover heater and sensor to avoid damage during using
Once-forming stainless steel lid
Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
Independent temperature-limiting alarm system
temperature error alarm
Test tube holder can be placed

SPECIFICATIONS

Model	BBSW-3005
Capacity	33 L
Temperature range	10-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1500 W
Interior Dimension	440Wx300Dx250H mm
Exterior Dimension	710Wx410Dx710H mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com