



## SHAKING WATER BATH BBSW-3005

## SHAKING WATER BATH BBSW-3005

### BBSW-3005 SHAKING WATER BATH



- A stamping molding stainless steel tank, easy to clean
- LCD screen, multiple data display with timing function, easy to operate
- Stainless steel shelves cover heater and sensor to avoid damage during using
- Once-forming stainless steel lid
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
- Independent temperature-limiting alarm system
- temperature error alarm
- Test tube holder can be placed

## SPECIFICATIONS

Model	BBSW-3005
Capacity	33 L
Temperature range	10-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1500 W
Interior Dimension	440Wx300Dx250H mm
Exterior Dimension	710Wx410Dx710H mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)