



## SHAKING WATER BATH BBSW-4001

# SHAKING WATER BATH BBSW-4001

## BBSW-4001 SHAKING WATER BATH



Unique and novel appearance, small size, powerful functions. All surface that contact with water adopt stainless steel material and PC material, with strong corrosion resistance

Using a once punch-formed stainless steel inner liner, with unique sealing structure to eradicate water leakage; Three groups of water tanks postposition independent drain valves separately, can easily emptying the sewage after clean the inner container

Using microcomputer program PID automatic temperature control, uniform temperature control, high precision

Built-in water level sensor, to prevent the hidden dangers of dry burning, safe and reliable

Dual digital tube digital holographic display temperature and timing time, parameter display is accurate and intuitionistic, with beep alarm prompt function

## SPECIFICATIONS

Model	BBSW-4001
Temperature Range	RT.+ 5 °C ~ 100 °C
Temp. accuracy	± 0.5 °C
Temp. Uniformity @ 37°C	± 0.5 °C
Number of holes	3
Workspace	276 × 118 × 98 mm × 1
Volume	15 L × 1
Dimension	W.570 × D.365 × H.465 mm
Net weight	13.6 kg
Time range	1 ~ 99 h 59 min (0 is forever)
Voltage	AC 220 V 50/60 Hz
Power	2200 W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)