



## WATER BATH BBWA-3101

## WATER BATH BBWA-3101

### BBWA-3101 WATER BATH

A stamping molding stainless steel tank, easy to clean  
LCD screen, multiple data display with timing function, easy to operate  
Stainless steel shelves cover heater and sensor to avoid damage during using  
Once-forming stainless steel lid  
Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time  
Independent temperature-limiting alarm system  
temperature error alarm  
Test tube holder can be placed

## SPECIFICATIONS

Model	BBWA-3101
Capacity	2 holes $\Phi 112$ mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	500 W
Temperature Stability	$\pm 0.3^{\circ}\text{C}$
Temp Alarm	$\pm 2^{\circ}\text{C}$
Interior dimension	130Wx280Dx150H mm
Exterior dimension	396Wx250Dx260H mm
Timing Range	1~5999 min



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)