



WATER BATH BBWA-3102

WATER BATH BBWA-3102

BBWA-3102 WATER BATH



- A stamping molding stainless steel tank, easy to clean
- LCD screen, multiple data display with timing function, easy to operate
- Stainless steel shelves cover heater and sensor to avoid damage during using
- Once-forming stainless steel lid
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
- Independent temperature-limiting alarm system
- temperature error alarm
- Test tube holder can be placed

SPECIFICATIONS

Model	BBWA-3102
Capacity	4 holes Φ 92 mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	1000 W
Temperature Stability	$\pm 0.3^{\circ}\text{C}$
Temp Alarm	$\pm 2^{\circ}\text{C}$
Interior dimension	220Wx280Dx150H mm
Exterior dimension	396Wx330Dx260H mm
Timing Range	1~5999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com