



WATER BATH BBWA-3102

WATER BATH BBWA-3102

BBWA-3102 WATER BATH



A stamping molding stainless steel tank, easy to clean
LCD screen, multiple data display with timing function, easy to operate
Stainless steel shelves cover heater and sensor to avoid damage during using
Once-forming stainless steel lid
Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
Independent temperature-limiting alarm system
temperature error alarm
Test tube holder can be placed

SPECIFICATIONS

Model	BBWA-3102
Capacity	4 holes $\Phi 92$ mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	1000 W
Temperature Stability	$\pm 0.3^{\circ}\text{C}$
Temp Alarm	$\pm 2^{\circ}\text{C}$
Interior dimension	220Wx280Dx150H mm
Exterior dimension	396Wx330Dx260H mm
Timing Range	1~5999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com