



WATER BATH BBWA-3105

WATER BATH BBWA-3105

BBWA-3105 WATER BATH

A stamping molding stainless steel tank, easy to clean
LCD screen, multiple data display with timing function, easy to operate
Stainless steel shelves cover heater and sensor to avoid damage during using
Once-forming stainless steel lid
Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
Independent temperature-limiting alarm system
temperature error alarm
Test tube holder can be placed

SPECIFICATIONS

| | |
|------------------------|-------------------|
| Model | BBWA-3105 |
| Capacity | 2holes+4holes |
| Temperature range | RT+5-100°C |
| Electrical Requirement | AC 220 50 Hz |
| Power Consumption | 500W+1000W |
| Temperature Stability | ±0.3°C |
| Temp Alarm | ±2°C |
| Interior dimension | 130Wx280Dx150H mm |
| Exterior dimension | 526Wx395Dx260H mm |
| Timing Range | 1~5999 min |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com