



WATER BATH BBWA-3107

WATER BATH BBWA-3107

BBWA-3107 WATER BATH



A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed

SPECIFICATIONS

Model	BBWA-3107
Capacity	11 L
Temperature range	RT+5-80°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	800 W
Temperature Stability	±0.2°C
Temp Alarm	0.1°C
Interior dimension	300Wx240Dx200H mm
Exterior dimension	480Wx300Dx480H mm
Timing Range	1~5999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com