



WATER BATH BBWA-3203

WATER BATH BBWA-3203

BBWA-3203 WATER BATH



Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature

SPECIFICATIONS

Model	BBWA-3203
Capacity	12 L
Electrical Requirement	220 V 50 HZ
Power Consumption	400 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	300Wx240Dx160H mm
Exterior dimension	460Wx280Dx190H mm
Timing Range	1~999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com