



WATER BATH BBWA-3204

WATER BATH BBWA-3204

BBWA-3204 WATER BATH



Microprocessor control with timing function
Digital display
Audible and visible alarm for over temperature

SPECIFICATIONS

Model	BBWA-3204
Capacity	22 L
Electrical Requirement	220 V 50 HZ
Power Consumption	600 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	450Wx300Dx190H mm
Exterior dimension	610Wx340Dx260H mm
Timing Range	1~999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com