



## WATER BATH BBWA-3205

# WATER BATH BBWA-3205

## BBWA-3205 WATER BATH



Microprocessor control with timing function  
Digital display  
Audible and visible alarm for over temperature

## SPECIFICATIONS

Model	BBWA-3205
Capacity	30 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	600Wx300Dx190H mm
Exterior dimension	760Wx340Dx280H mm
Timing Range	1~999 min



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)