



WATER BATH BBWA-3206

WATER BATH BBWA-3206

BBWA-3206 WATER BATH



- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature

SPECIFICATIONS

Model	BBWA-3206
Capacity	2.1 L x 3
Electrical Requirement	220 V 50 HZ
Power Consumption	750 W
Temperature Range	RT+5-99°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	150Wx125Dx110H mm
Exterior dimension	490Wx245Dx310H mm
Timing Range	1~999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com