



BENCHTOP HIGH SPEED CENTRIFUGE BCBH-302


BENCHTOP HIGH SPEED CENTRIFUGE BCBH-302

Designed around your applications to provide you competent sample processing and reliable results. Microcomputer programmable with excellent temperature controls and low noise operation maximize your productivity. Highly efficient with low maintenance requirements, it is an ideal separation tool for multiple research applications.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

BCBH-302 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running

Automatically electric lid lock, super speed, over temperature protection and imbalance protection

The centrifuge body is made of high quality steel, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable

Microprocessor control, TD5A with digital display and TD5C with LCD display which indicate the speed, time, RCF in operation, speed raising and reducing quick, operate simply

10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice

Three layers protection of steel structure and get the ideal centrifugation result

SPECIFICATIONS

Model	BCBH-302
Maximum Capacity (No of tubes x Vol.)	6x100 ml
Maximum Speed	21000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	28030xg
Time Range	1~99 min
Overall Dimension	513x370x320 mm
Noise Level	<60 dB(A)
Net. Weight	42 kg
Power Supply	AC 220V 50HZ 10A
Package	carton box

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301008006	Angle rotor	12x1.5/2 ml	21000	28030xg
2301008007	Angle rotor	24x1.5/2 ml	17000	26460xg
2301008008	Angle rotor	16x5 ml	16000	22020xg
2301008009	Angle rotor	6x30 ml	14000	19060xg
2301008010	Angle rotor	12x10 ml	15000	22680xg
2301008011	Angle rotor	12x15 ml	10000	11840xg

2301008012	Angle rotor	6x10 ml	16000	21500xg
2301008013	Angle rotor	6x50 ml	12000	16050xg
2301008014	Angle rotor	8x20 ml	15000	22680xg
2301008015	Angle rotor	6x70 ml	12000	15570xg
2301008016	Angle rotor	6x100 ml	10000	11380xg
2300908024	Bucket rotor	2x3x48 well	4000	2300xg
2301008017	Angle rotor	40x0.5 ml	15000	22920xg
2301008018	Angle rotor	30x1.5/2 ml	13500	19340xg
2301008019	Angle rotor	12x7 ml	16000	21380xg
2301008020	Angle rotor	4x100 ml	12000	14850xg
2301008021	Angle rotor	30x0.5 ml	18000	26660xg
2301008022	Angle rotor	48x1.5/2 ml	13000	17930xg
2301008023	Vertical rotor	16x5 ml	14000	12660xg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com