



BENCHTOP HIGH SPEED REFRIGERATED CENTRIFUGE BCBHR-304

BENCHTOP HIGH SPEED REFRIGERATED CENTRIFUGE

BCBHR-304

Engineered to maximize productivity and reliably deliver a wide range of high speed separations. Stainless steel interiors, integrated refrigeration system and advanced temperature controls make it flexible for your evolving applications. Safety features and low noise operation make it an ideal partner for your day to day separations. Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Floor Type Refrigerated Centrifuge, Benchtop Centrifuge, Laboratory Tabletop Centrifuge.

BCBHR-304 BENCHTOP HIGH SPEED REFRIGERATED CENTRIFUGE



- Microcomputer controlled programmable operation with push-button operation
- TFT true-colour LCD wide screen displays set and run conditions for easy monitoring
- Triple protective layers of steel bushings
- Brushless DC Motor ensures low maintenance, high efficiency and constant speed of rotation
- Fully stainless steel structure and a stainless steel cavity makes the machine safe and more efficient
- Acceleration/Deceleration takes around 30 sec to 1 min
- Multiple layers of shock absorbers with auto balancing function
- System is complemented by multiple new accessories, including tube racks, centrifuge adapters and rotors
- Safety features: Self diagnostic, Rotor unbalance detector, Lid protection, Over speed protection, Over temperature protection

SPECIFICATIONS

Model	BCBHR-304
Maximum Capacity (No of tubes x Vol.)	4x100 ml
Temperature Range	-20°C to +40°C
Temperature Accuracy	±1°C
Maximum Speed	16800rpm
Speed Increment	±50 rpm
Maximum RCF	20320×g
Time Range	0-99 h 59 min
Overall Dimension	710x630x350 mm
Program	100 user programs
Control	Microcomputer
Display	LCD
Noise Level	≤65dB
Weight (Net/Gross)	75/95 kg
Power	450 W
Power Supply	220V, 60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301314006	Fixed Angle Rotor	12x1.5/2.0 ml	16800	19600 _{xg}
2301314007	Fixed Angle Rotor	24x1.5/2.0 ml	14000	18180 _{xg}
2301314008	Fixed Angle Rotor	18x1.5 ml	14000	13600 _{xg}
2301314009	Fixed Angle Rotor	30x1.5 ml	11000	12900 _{xg}
2301314010	Fixed Angle Rotor	48x0.5 ml	12000	13800 _{xg}
2301314011	Fixed Angle Rotor	10x 5 ml	13000	11700 _{xg}
2301314012	Fixed Angle Rotor	12x10ml	12000	13400 _{xg}
2301314013	Fixed Angle Rotor	12x5ml	15000	15600 _{xg}
2301314014	Fixed Angle Rotor	8x15 ml	12000	14800 _{xg}
2301314015	Fixed Angle Rotor	6x50 ml	11000	11600 _{xg}
2301314016	Fixed Angle Rotor	4x100 ml	9000	8500 _{xg}
2301314017	Swinging-bucket Rotor	4x100 ml	4000	2100 _{xg}
2301314018	Swinging-bucket Rotor	12x10 ml	4000	2330 _{xg}
2301314019	Swinging-bucket Rotor	4x50 ml	4000	2100 _{xg}
2301314020	Swinging-bucket Rotor	4x10 ml	6000	3200 _{xg}
2301314021	Swing Micro-Plate Rotor	2x2x48 2x2x48 places	3000	1400 _{xg}



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com