



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-104

## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-104

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-104 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



## SPECIFICATIONS

Model	BCBT-104
Capacity	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	1090x460x1120 mm
Package Size	1210x680x1285 mm
Gross Weight	117 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)