



## STRONG ACID AND ALKALI STORAGE CABINET BCBT-201

## STRONG ACID AND ALKALI STORAGE CABINET BCBT-201

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-201 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable

## SPECIFICATIONS

Model	BCBT-201
Capacity	4 gallon (15L)
Material	pp
Adjustable Shelf	Optional
Door	Single / Manual
External Size	430x430x500 mm
Package Size	515x515x685 mm
Gross Weight	16 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)