



STRONG ACID AND ALKALI STORAGE CABINET BCBT-203

STRONG ACID AND ALKALI STORAGE CABINET BCBT-203

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-203 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable

SPECIFICATIONS

Model	BCBT-203
Capacity	20 gallon (75L)
Material	pp
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	910x590x890 mm
Package Size	995x675x1035 mm
Gross Weight	46 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com