



# STRONG ACID AND ALKALI STORAGE CABINET BCBT-204

## STRONG ACID AND ALKALI STORAGE CABINET BCBT-204

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-204 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable

## SPECIFICATIONS

Model	BCBT-204
Capacity	28 gallon (105L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	910x600x900 mm
Package Size	995x685x1035 mm
Gross Weight	48 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)