



STRONG ACID AND ALKALI STORAGE CABINET BCBT-205

STRONG ACID AND ALKALI STORAGE CABINET BCBT-205

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-205 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable

SPECIFICATIONS

Model	BCBT-205
Capacity	30 gallon (114L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1120mm
Package Size	1175x545x1255 mm
Gross Weight	58 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com