



## SLIDE CABINET BCBT-301-A

## SLIDE CABINET BCBT-301-A

Biolab's Slide Cabinet is equipped with sliding rails underneath the drawers which offers smooth sliding while inside hidden locks prevent unexpected opening. The interior space of each drawer is divided into two sections to make filing more convenient.

Used in Pathological laboratories, Research institutions, Universities, Cytology and Histology applications.

### BCBT-301-A SLIDE CABINET



Sliding rails underneath the drawers offer smooth sliding while inside hidden locks prevent unexpected opening.

The interior space of each drawer is divided into 6 sections to make filing more convenient.

## SPECIFICATIONS

Model	BCBT-301-A
Drawers	72 drawers (12 layers and 6 drawers in each layer) (each can be individually locked)
Capacity	64,800 microscope slides with 900 in each drawer
Base Material	SPCC1.2 cold-roll steel sheets
Main body and Drawers Material	SPCC 0.8 cold-roll steel sheets
Dimensions	450x478x1620 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)