



CLIMATIC CHAMBER BCCL-203

CLIMATIC CHAMBER BCCL-203

Our chambers provide excellent temperature and lighting control for regulated growth of plants and other products. Its responsive humidification sensor and energy efficient design make it a preferred choice for a wide range of applications.

Used in Supports plant and insect growth, tissue cultivation, antibiotic cultivation and preservation and other controlled temperature experiments..

Also known as Humidity Chamber, Laboratory Climatic Chamber.

BCCL-203 CLIMATIC CHAMBER



- Forced convection mode
- 30 Stages programmable Microcomputer PID
- Imported environmental protection insulation materials
- Excellent Imported temperature and humidity Sensor
- SUS304 stainless steel interior material
- ALLCOLDTM: Refrigeration Technology
- ALLIGENTTM: Humidification System
- ALLFLOWTM: Perfect Air Current Cycling

SPECIFICATIONS

Model	BCCL-203
Capacity	400 L
Temperature Range	10°C~65°C(Light on)/0°C~60°C(Light off)
Temperature Uniformity	±1.5 (10-40°C)
Temperature Fluctuation	±1°C
Temperature Resolution	0.1°C
Using Environment Temperature	Ambient temp.10~30°C
Using Environment Humidity	<70%
Illumination Range	0-15000 (5 stages adjustable)
Timer Range	0-99h/0-9999mins
Shelves	4
Inner Dimension	1148 x 554 x 610 mm
Overall Dimension	1828 x 783 x 905 mm
Weight (Net/Gross)	194 kg
Power	1550 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com