



## WATER CHILLER BCHI-410

# WATER CHILLER BCHI-410

## BCHI-410 WATER CHILLER



- P.I.D temperature controller provides accurate and reliable temperature control
- Large LCD display screen and interface provides for user-friendly operation
- Preset On/Off function
- Independent circulating pump switch for easy starting/stopping the cycle
- Easy to set Adjustable Timer. (1 minute to 5,999 minutes)
- Maintenance-free operation with easy to clean surface
- Both heater and bath chamber are made of corrosion-resistant stainless steel
- No angle in bath chamber, easy to clean
- Uniquely designed heating method can up to 80 °C

## SPECIFICATIONS

Model	BCHI-410
Storage tank maximum capacity	20 L
Temperature range	-20-40°C
Ambient temperature range	+5-30°C
Temperature accuracy	±3°C
Cooling capacity at 20°C	3.52 Kw
Cooling capacity at 0°C	2.96 Kw
Cooling capacity at -20°C	1.3 Kw
Refrigerant	R404A
Security features	Delay, leakage, overcurrent, overvoltage
Total power	3800 W
Heating power	2000 W
Power requirements	AC 220 V ± 10% /50 HZ
Pump flow max.	17 L/min
Maximum head	2.5 m
Inlet/Outlet pipe diameter	Φ20
Noise level	≤55
Liquid tank opening / depth	Φ300x300 mm
Dimensions	430Wx600Dx840H mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)