



## HIGH TEMPERATURE CIRCULATOR BCHT-101

## HIGH TEMPERATURE CIRCULATOR BCHT-101

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water bath.

## BCHT-101 HIGH TEMPERATURE CIRCULATOR



Microprocessor temperature controller

Easy maintenance

LCD display with touch key panel for temperature setting

Over temperature alarm system

## SPECIFICATIONS

Model	BCHT-101
Capacity	5 L
Temperature Range	RT+5°C to 95°C
Temperature Accuracy	±0.2°C
Pump Flow Rate	6 L/min
Top Opening	130x130 mm
Inner Dimension (WxDxH)	240x140x150 mm
Outer Dimension (WxDxH)	390x350x520 mm
Weight	9 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)