



HIGH TEMPERATURE CIRCULATOR BCHT-102

HIGH TEMPERATURE CIRCULATOR BCHT-102

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water bath.

BCHT-102 HIGH TEMPERATURE CIRCULATOR



- Microprocessor temperature controller
- Easy maintenance
- LCD display with touch key panel for temperature setting
- Over temperature alarm system

SPECIFICATIONS

Model	BCHT-102
Capacity	15 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	235x160 mm
Inner Dimension (WxDxH)	300x240x200 mm
Outer Dimension (WxDxH)	430x380x570 mm
Weight	12.84 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com