



HIGH TEMPERATURE CIRCULATOR BCHT-103

HIGH TEMPERATURE CIRCULATOR BCHT-103

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water bath.

BCHT-103 HIGH TEMPERATURE CIRCULATOR



Microprocessor temperature controller

Easy maintenance

LCD display with touch key panel for temperature setting

Over temperature alarm system

SPECIFICATIONS

| | |
|-------------------------|-----------------|
| Model | BCHT-103 |
| Capacity | 25 L |
| Temperature Range | RT+5°C to 100°C |
| Temperature Accuracy | ±0.05°C |
| Pump Flow Rate | 6 L/min |
| Top Opening | 235x160 mm |
| Inner Dimension (WxDxH) | 280x250x300 mm |
| Outer Dimension (WxDxH) | 460x440x670 mm |
| Weight | 20 kg |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com