



HIGH TEMPERATURE CIRCULATOR BCHT-103

HIGH TEMPERATURE CIRCULATOR BCHT-103

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water bath.

BCHT-103 HIGH TEMPERATURE CIRCULATOR



- Microprocessor temperature controller
- Easy maintenance
- LCD display with touch key panel for temperature setting
- Over temperature alarm system

SPECIFICATIONS

Model	BCHT-103
Capacity	25 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	235x160 mm
Inner Dimension (WxDxH)	280x250x300 mm
Outer Dimension (WxDxH)	460x440x670 mm
Weight	20 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com