



LOW TEMPERATURE CIRCULATOR BCLT-2303

LOW TEMPERATURE CIRCULATOR BCLT-2303

BCLT-2303 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller with LCD screen ensures precise and reliable control, easy to operate

No angle in bath chamber, easy to clean

With interface to external water bath

Easy to set adjustable timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel

Drain valve makes emptying of bath water fast and easy for cleaning and moving

Programming setting function with 7 periods and 9 steps for each period, which means there are 63 programmable steps in total

SPECIFICATIONS

Model	BCLT-2303
Capacity	22.5 L
Temperature range	RT+5 - 150°C
Precision	±0.1
Liquid tank opening depth	330x300/150 mm
Power requirements	AC 220 V 50 HZ
Pump	8 L/min
Power Consumption	1050 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com