



LOW TEMPERATURE CIRCULATOR BCLT-2801

LOW TEMPERATURE CIRCULATOR BCLT-2801

BCLT-2801 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller
Audible and visible alarm for temperature and water level
R134a refrigerant, imported compressor
With interface to external water bath
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

SPECIFICATIONS

Model	BCLT-2801
Capacity	6.7 L
Temperature range	RT+5-100°C
Precision	±0.1
Interior Dimension	150x160x150
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com