



LOW TEMPERATURE CIRCULATOR BCLT-2804

LOW TEMPERATURE CIRCULATOR BCLT-2804

BCLT-2804 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller
Audible and visible alarm for temperature and water level
R134a refrigerant, imported compressor
With interface to external water bath
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

SPECIFICATIONS

Model	BCLT-2804
Capacity	14.5 L
Temperature range	RT+5-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com