



## LOW TEMPERATURE CIRCULATOR BCLT-2805

# LOW TEMPERATURE CIRCULATOR BCLT-2805

## BCLT-2805 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller  
Audible and visible alarm for temperature and water level  
R134a refrigerant, imported compressor  
With interface to external water bath  
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

## SPECIFICATIONS

|                         |             |
|-------------------------|-------------|
| Model                   | BCLT-2805   |
| Capacity                | 14.5 L      |
| Temperature range       | RT+5-100°C  |
| Precision               | ±0.2        |
| Interior Dimension      | 240x170x200 |
| Electrical Requirements | 220 V 50 Hz |
| Pump (flux)             | 8 L/min     |
| Power Consumption       | 1050 W      |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)