



## LOW TEMPERATURE CIRCULATOR BCLT-2806

## LOW TEMPERATURE CIRCULATOR BCLT-2806

### BCLT-2806 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller  
Audible and visible alarm for temperature and water level  
R134a refrigerant, imported compressor  
With interface to external water bath  
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

## SPECIFICATIONS

Model	BCLT-2806
Capacity	4.5 L
Temperature range	-10-100°C
Precision	±0.2
Interior Dimension	150x160x150
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	2300 W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)