



LOW TEMPERATURE CIRCULATOR BCLT-2807

LOW TEMPERATURE CIRCULATOR BCLT-2807

BCLT-2807 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller
Audible and visible alarm for temperature and water level
R134a refrigerant, imported compressor
With interface to external water bath
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

SPECIFICATIONS

| | |
|-------------------------|-------------|
| Model | BCLT-2807 |
| Capacity | 4.5 L |
| Temperature range | -20-100°C |
| Precision | ±0.2 |
| Interior Dimension | 150x160x150 |
| Electrical Requirements | 220 V 50 Hz |
| Pump (flux) | 8 L/min |
| Power Consumption | 2300 W |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com