



## LOW TEMPERATURE CIRCULATOR BCLT-2811

# LOW TEMPERATURE CIRCULATOR BCLT-2811

## BCLT-2811 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller  
Audible and visible alarm for temperature and water level  
R134a refrigerant, imported compressor  
With interface to external water bath  
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

## SPECIFICATIONS

Model	BCLT-2811
Capacity	13 L
Temperature range	-10-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	2300 W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)