



LOW TEMPERATURE CIRCULATOR BCLT-2813

LOW TEMPERATURE CIRCULATOR BCLT-2813

BCLT-2813 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller
Audible and visible alarm for temperature and water level
R134a refrigerant, imported compressor
With interface to external water bath
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

SPECIFICATIONS

Model	BCLT-2813
Capacity	13 L
Temperature range	-40-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	3100 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com