



## LOW TEMPERATURE CIRCULATOR BCLT-2815

# LOW TEMPERATURE CIRCULATOR BCLT-2815

## BCLT-2815 LOW TEMPERATURE CIRCULATOR



Microprocessor temperature controller  
Audible and visible alarm for temperature and water level  
R134a refrigerant, imported compressor  
With interface to external water bath  
RS485 connector is option which can connect computer to record the parameters and the variations of temperature

## SPECIFICATIONS

Model	BCLT-2815
Capacity	≤50 L
Temperature range	-100-200°C
Precision	0.1
Interior Dimension	130x150x330
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)