



## LOW TEMPERATURE CIRCULATOR BCLT-2903

## LOW TEMPERATURE CIRCULATOR BCLT-2903

### BCLT-2903 LOW TEMPERATURE CIRCULATOR



High precision temperature controller

## SPECIFICATIONS

Model	BCLT-2903
Capacity	13 L
Temperature range	-40-100°C
Precision	±0.02
Peristome Dimension	240x170x200
Pump(flux)	15 L/min
Power	2850 W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)