



MINI CENTRIFUGE BCMI-101

MINI CENTRIFUGE BCMI-101

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

BCMI-101 MINI CENTRIFUGE

Compact, corrosion resistant design provides low maintenance with longer life

Micro DC motor and carbon-fibre rotor combined with silence technology reduces the operating noise

Safety features: Over speed and over temperature protection

SPECIFICATIONS

Model	BCMI-101
Maximum Capacity (No of tubes x Vol.)	6x1.5/0.5/0.2 ml, 0.2 mlx2x8
Maximum Speed	4000 rpm
Speed Increment	±50 rpm
Maximum RCF	900xg
Time Range	None
Overall Dimension	189x161x122 mm
Display	LED
Noise Level	< 48 dB(A)
Power Supply	220 V, 60 Hz

ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2302306006	Rotor	1.5/0.5/0.2 mlx6	4000	900xg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com