



MINI CENTRIFUGE BCMI-101

MINI CENTRIFUGE BCMI-101

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

BCMI-101 MINI CENTRIFUGE



Compact, corrosion resistant design provides low maintenance with longer life
Micro DC motor and carbon-fibre rotor combined with silence technology reduces the operating noise
Safety features: Over speed and over temperature protection

SPECIFICATIONS

| | |
|---------------------------------------|-----------------------------|
| Model | BCMI-101 |
| Maximum Capacity (No of tubes x Vol.) | 6x1.5/0.5/0.2 ml, 0.2mlx2x8 |
| Maximum Speed | 4000 rpm |
| Speed Increment | ±50 rpm |
| Maximum RCF | 900xg |
| Time Range | None |
| Overall Dimension | 189x161x122 mm |
| Display | LED |
| Noise Level | < 48 dB(A) |
| Power Supply | 220 V, 60 Hz |

ACCESSORIES

| Accessory Code | Name | Description | RPM | RCF _{xg} |
|----------------|-------|-----------------|------|-------------------|
| 2302306006 | Rotor | 1.5/0.5/0.2mlx6 | 4000 | 900xg |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com