



MINI CENTRIFUGE BCMI-302

MINI CENTRIFUGE BCMI-302

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

BCMI-302 MINI CENTRIFUGE

Micro brushless DC motor provides quiet, maintenance-free operation

Compact design with small footprint

Optimized air flow minimizes sample warming and protects temperature-sensitive samples

UV and flame-retardant composite material ensures high quality

Triple flip switch function opens the lid automatically at the end of the run to prevent sample warming and allow easy access to samples

SPECIFICATIONS

Model	BCMI-302
Maximum Capacity (No of tubes x Vol.)	6x1.5ml/0.5ml/0.2ml
Maximum Speed	12000 rpm
Speed Increment	1000 rpm
Maximum RCF	7450xg
Time Range	1 - 99 min
Overall Dimension	220x210x120 mm
Packing Size	153x128x104 mm
Display	LED
Noise Level	< 56 dB(A)
Weight (Net/Gross)	5 kg
Power	10 W
Power Supply	85~245 V/ 60 Hz

ACCESSORIES

Accessory Code	Name	Description
2302507006	Rotor	8x1.5 ml

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
2302507007	Rotor	4x8x0.2 ml
2302507008	Rotor - Adaptor	8x0.5 ml
2302507009	Rotor - Adaptor	8x0.2 ml



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com