



## MINI CENTRIFUGE BCMI-404

## MINI CENTRIFUGE BCMI-404

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

## BCMI-404 MINI CENTRIFUGE



Compact design takes up less than six inches of bench space

Easy open clear colored lids

Simply close the lid for acceleration and open the lid for deceleration, or alternatively use the on/off switch to control

Ideal for quick spin downs of micro tubes and PCR tubes

Safety Features: Over-speed detection, Over-temperature detection

## SPECIFICATIONS

Model	BCMI-404
Maximum Capacity (No of tubes x Vol.)	2 mlx24, 0.5mlx36, PCR 8 serial tubes
Maximum Speed	200-15000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	21380xg
Time Range	30 s- 99 min
The shortest Acceleration time	25 sec
The shortest Deceleration time	25 sec
Overall Dimension	280x364x266 mm
Display	LCD
Noise Level	≤ 58 dB(A)
Weight (Net/Gross)	12 kg
Power Supply	single phase, AC220-240V/50Hz/60Hz,5A; AC110V-120V,50hz60hz,8A

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
2302609006	Aluminum alloy rotor kit	Rotor kit with lid; fixing clips and O'ring; 15000rpm; 5mlx18
2302609007	Aluminum alloy rotor kit	Rotor kit with lid; fixing clips and O'ring; 15000rpm; 5mlx12 (eppendorf tube)
2302609008	Aluminum alloy rotor kit	Rotor kit with lid; fixing clips and O'ring; 15000rpm; 2ml/1.5mlx24

Accessory Code	Name	Description
2302609010	Aluminum alloy rotor kit	Rotor kit with lid; fixing clips and O'ring; 15000rpm; 4xPCR 8
2302607006	Rotor adapter	2ml to 0.2ml rotor adapter; 24pcs/pk
2302607007	Rotor adapter	2ml to 0.5ml rotor adapter; 24pcs/pk

2302609009	Aluminum alloy rotor kit	Rotor kit with lid; fixing clips and O'ring; 15000rpm; 0.5mlx36
------------	--------------------------	-----------------------------------------------------------------

2302609011	Rotor adapter	0.5ml to 0.2ml rotor adapter; 36pcs/pk
------------	---------------	----------------------------------------



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)