



## MINI CENTRIFUGE BCMI-502

## MINI CENTRIFUGE BCMI-502

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

## BCMI-502 MINI CENTRIFUGE



Large radius rotor with 8 tube holes, centrifugal force is 1.5 times bigger than ordinary 6-hole rotor at the same rolling speed

The cover press-button is outward designed with start/stop function. Press the button, the cover will open automatically to 95°, which is very convenient for one-hand operation

The rotor is clip-on fixed. Very convenient for rotors exchanging yet fixed solid reliable

LED displays speed (MINI-10K+) and timing

Working is silent and stable. Free-maintenance motor, long-life and safe

Smart appearance. Multi-purpose to meet different experimental demands

## SPECIFICATIONS

Model	BCMI-502
Rotor Capacity (No of tubes x Vol.)	8x1.5/2.0 ml, 8x0.5 ml, 16x0.2 ml(single tube or 2x8x0.2 ml strips). (With 0.5 ml and 0.2 ml tube adapter)
Speed Range	4000 rpm
RCF	Approx.700 g
Time Range	1s ~9999 min or continuous
Noise Level	≤45 dB
Voltage	AC 220 V / AC 110 V, 50/60 Hz
Power	≤ 45 W
Fuse	250 V, 1 A, Φ5x20
Overall Dimension	W.156 x D.176 x H.121 mm
Net Weight	1.5 kg

## ACCESSORIES

Accessory Code	Name
2302707006	Rotor
2302707007	Strip Rotor
2302707008	Tube adaptors



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)