



## OIL BATH CIRCULATOR BCOB-104

## OIL BATH CIRCULATOR BCOB-104

Designed with intelligent control system and wide temperature range, oil bath circulator is an ideal for industrial and laboratory applications. Available with varying capacities and circulating pump requirements, it has a user friendly interface and an optional water circulation cooling device for controlling exothermic reactions.

Used in High temperature organic distillation, External circulation, Calibration tests, High temperature experiments, Laboratory.

Also known as Cooling Bath Circulator, Refrigerated Circulator, Laboratory Refrigerated Bath Circulator.

## BCOB-104 OIL BATH CIRCULATOR

Advanced inner circulation and outer circulation pump

LED dual-window digital display

High effective thermal insulation materials

## SPECIFICATIONS

Model	BCOB-104
Capacity	50 L
Temperature Range	RT +5 to 300°C
Temperature Accuracy	±0.1°C
Pump Flow Rate	8 L/min
Inner Dimension (WxDxH)	400x350x350 mm
Control	Microprocessor Controller
Display	LED
Power Supply	220 V



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)