



## FREEZE DRYER BFFT-303-A

## FREEZE DRYER BFFT-303-A

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## BFFT-303-A FREEZE DRYER



## SPECIFICATIONS

Model	BFFT-303-A
Type	Manifold Type
Condenser Temperature	-110°C
Freeze Drying area	0.12 m <sup>2</sup>
Bulk capacity	1.2 L
Tray	φ200mm, 4pcs
Round bottom flask	1000ml,500ml ,250ml



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)