



FREEZE DRYER BFFT-317-B

FREEZE DRYER BFFT-317-B

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

BFFT-317-B FREEZE DRYER



SPECIFICATIONS

Model	BFFT-317-B
Type	Manifold type Electric heating shelf
Condenser Temperature	-110°C
Freeze Drying area	0.36 m ²
Bulk capacity	3.6 L
Tray	φ240mm, 8pcs
Round bottom flask	1000ml,500ml ,250ml ,100ml,50ml



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com