



## UPRIGHT FREEZER DOUBLE DOOR BFRZ-25-403

## UPRIGHT FREEZER DOUBLE DOOR BFRZ-25-403

High performance refrigerators featuring microprocessor control for temperature accuracy and ergonomic design for easy usage. Designed for space saving storage, it provides outstanding thermal performance along with safety features.

Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## BFRZ-25-403 UPRIGHT FREEZER DOUBLE DOOR



## SPECIFICATIONS

Model	BFRZ-25-403
Capacity	305 L
Temperature Range	-10°C to 25°C
Alarm	High&low temp, control error, power failure
External Dimension	1135x725x890 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)